



Unique Reverse Phase Interaction FluoroSep-RP Phenyl / HS

# Application News

701 South Route 73 • West Berlin, New Jersey 08091-2621 • (856) 753-8400 • Fax (856) 753-8484

## Separation of Amino Phenols

Column: FluoroSep-RP Phenyl / HS 5 $\mu$  120Å 25cm X 4.6mm  
Catalog Number: 155211-FSP/HS

### Test Conditions

Mobile Phase: 88% Sodium Citrate 50mM pH=6.4 / 15% ACN  
Flow Rate: 1.0 mL/min  
Detection: UV at 234nm  
Injection vol.: 5  $\mu$ L  
Backpressure: 160 bar  
Test compounds: Amino Phenols  
Sample: 125 $\mu$ g/mL per component

